


## After you receive cell lines

Transfer quickly from the package to storing place or culturing.

Frozen cells in above -65°C will be rapidly damaged.



# First Culturing Guidance

#### Media and equipment

- Use a 25cm<sup>2</sup> flask or a 6cm dish (or smaller equipment if it is specified in data sheet) to avoid excessive dilution.
- Prepare the appropriate media specified in the data sheet of the web catalog.

If you have any question, please inquire with JCRB Cell Bank before you start culturing.

[Q&A: the reference page in JCRB web site] http://cellbank.nibiohn.go.jp/english/other\_e/other\_faq\_e.html

#### Cell lines should be thawed rapidly

[Handle one by one. Don't thaw two or more ampules together.]

- 1. Take out one ampule from the package. (Use safety gloves and a face shield)
- 2. Immediately put the ampule in warm (not higher than 37°C) water, and thaw the content within 2 min by shaking.

### Seed at 25 cm<sup>2</sup> flask or 6 cm dish

[Under aseptic conditions]

- 3. Prepare 10 ml of the specified medium in a centrifuge tube.
- 4. Sterilize the ampule by permeating gauze with 70% ethanol or cationic detergent sterilizer
- 5. Wrap the ampule with sterilized gauze, and snap off the neck of the ampule with care.
- 6. Transfer the cell suspension to a centrifuge tube.
- 7. Centrifuge the mixture at 1000 rpm for 5 min, and discard the supernatant.
- 8. Without washing, re-suspend the cells in 5ml of the same medium. If especially the media volume is recommended in the data sheet, follow the instruction.
- Culture them in a 25 cm<sup>2</sup> flask or a 6 cm dish for tissue cultures. The cell density given in the attached data sheet is standard density for logarithmic growth phase. It's not the actual density on thawing.
- 10. Make sure of cell proliferation before proceeding to passage. *It may be needed for some days or a week for proliferation.* After growing well, make frozen stocks early.

